


multisub

HORIZONTAL PACKAGE DEALS



MINIPACK



MIDIPACK



CHOICEPACK



MAXIPACK



MULTISUB4PACK



SCREENPACK

ORDERING INFORMATION

- MINIPACK** **Mini Horizontal Electrophoresis Package**, includes: - BCMSMINIDUO - horizontal unit with 7x7cm and 7x10cm gel trays and combs (Page 6) and NANO300 - 300V, 400mA, 60W mini power supply (Page 49)
- MIDIPACK** **Midi Horizontal Electrophoresis Package**, includes: - BCMSMIDIDUO - horizontal unit with 10x7cm and 10x10cm gel trays and combs (Page 7) and NANO300 - 300V, 400mA, 60W mini power supply (Page 49)
- CHOICEPACK** **Choice Horizontal Electrophoresis Package**, includes: - BCMSCHOICETRIO - horizontal unit with 15x7cm, 15x10cm and 15x15cm gel trays and combs (Page 8) and MP-300 - 300V, 700mA, 150W midi power supply (Page 49)
- MAXIPACK** **Maxi Horizontal Electrophoresis Package** with extra gel tray, includes: - BCMSMAXIDUO - horizontal unit with 20x10cm and 20x20cm gel trays and combs (Page 9) and MP-300 - 300V, 700mA, 150W midi power supply (Page 49)
- SCREENPACK** **Screen Horizontal Electrophoresis Package**, includes: - BCMSSCREENTRIO - horizontal unit with 26x16cm, 26x24cm and 26x32cm gel trays and combs (Page 10) and MP-500 - 500V, 800mA, 300W maxi power supply (Page 49)
- MIDI96PACK** **Midi96 Horizontal Electrophoresis Package**, includes: - BCMSMIDI96 - horizontal unit with UV tray and 12x8-sample 1mm comb block (Page 12) and MP-300 - 300V, 700mA, 150W mini power supply (Page 49)
- MULTISUB4PACK** **MultiSUB-4 Horizontal Electrophoresis**, includes: - BCMULTISUB-4 - horizontal unit with four 8x12cm gel trays and combs (Page 13) and MP-300 - 300V, 700mA, 150W midi power supply (Page 49)