


# 4-gel tetrad systems

VERTICAL PACKAGES



TETRAD10



TETRAD10-CBS




TETRAD10-PS



TETRAD10-CBS-PS

## standard vertical electrophoresis



OMNI10-PS



OMNI10-CBS-PS



OMNI20-CBS-PS

### ORDERING INFORMATION

- TETRAD10** **OmniPAGE TETRAD Mini-Protein Electrophoresis Package** for 1mm self-cast gels. Includes: - BCVS10DSYS-CU – 10x10cm mini-vertical unit with casting base and external casting stand, plus 2x 1mm 12-sample combs, 2x plain glass plates with 1mm spacer and 2x notched plates (Page 32), VS10-12-1 – 2x additional 1mm 12-sample combs, VS10PGS1 – 1 pack of 2x plain glass plates with 1mm spacer, VS10NG – 1 pack of 2x notched glass plates with 1mm spacer, VS10NGS1 - 2 packs of 2x notched glass plates with 1mm spacer
- TETRAD10-CBS** **OmniPAGE TETRAD Mini-Protein Electrophoresis and Blotting Package** for 1mm self-cast gels and interchangeable 4-blot module, Includes: - TETRAD10 and VS10BI – 4-blot insert with 4 blotting cassettes and 8 fibre pads (Page 35)
- TETRAD10-PS** **OmniPAGE TETRAD** with standard Midi power supply option\*. Includes: - TETRAD10 - OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels and MP-300 - 300V, 700mA, 150W midi power supply (Page 53)
- TETRAD10-CBS-PS** **OmniPAGE TETRAD-CBS** with standard Midi power supply\*. Includes: - TETRAD10-CBS - OmniPAGE TETRAD10 Mini-Protein Electrophoresis Package for 1mm self-cast gels with blotting module and MP-300 - 300V, 700mA, 150W midi power supply (Page 53)
- OMNI10-PS** **OmniPAGE Mini Vertical** with standard Midi power supply option\*. Includes: - BCVS10DSYS - OmniPAGE Mini Vertical with casting base and MP-300 - 300V, 700mA, 150W midi power supply (Page 53)
- OMNI10-CBS-PS** **OmniPAGE Mini Vertical** with interchangeable 4-blot module and standard Midi power supply option\*. Includes: - OMNI10-PS - OmniPAGE Mini Vertical with casting base, VS10BI – 4-blot insert with 4 blotting cassettes and 8 fibre pads (Page 35) and MP-300V - 300V, 700mA, 150W midi power supply (Page 53)
- OMNI20-CBS-PS** **OmniPAGE WAVE Maxi Vertical\*** with interchangeable 3-blot module and standard Midi power supply option\*. Includes: - BCWAVEDYS OmniPAGE WAVE Maxi Vertical with casting base, VS20BI – 3-blot insert with 3 blotting cassettes and 6 fibre pads (Page 35) and MP-300 - 300V, 700mA, 150W midi power supply (Page 53)

\* 1) all packages available with High Voltage MP-500P Power Supply, 2) WAVE Maxi Tetrad systems also available, please contact for information and pricing